



Prime Comfort Active Hybrid Mattress

The Prime Comfort Active mattress combines the benefits of the Prime Comfort Plus high specification foam mattress with a dynamic surface of alternating air cells. This lets carers step patient care up or down based on changes in condition. The castellated foam layer provides enhanced pressure redistribution and keeps patients supported even when the cells are alternating or the mattress is deflated. Designed for community settings and patients at very high risk.

KEY PRODUCT FEATURES

- ✓ **Hybrid Foam & Air Design:** Castellated foam topper over 10 alternating cells for step up or step down care.
- ✓ **Firm 'U' Shaped Base:** High-density foam base and edges for roll prevention and to aid transfers.
- ✓ **Protected When Deflated:** Patients at very high risk stay protected even when the mattress is deflated.
- ✓ **Versatile Pump Modes:** Alternating, Constant Low Pressure, Max Inflate and Seated modes.

HIGHLIGHTS

- ★ **Hybrid Foam & Air Cell Design**
- ★ **Very High Risk with Pump**
- ★ **Comfort Control Settings**
- ★ **247kg Max User Weight**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Max Inflate with auto return to alternating for nursing procedures.
- Seated mode to support the patient when the backrest is elevated.
- Pump maintains set pressures to reduce carer intervention.
- No minimum user weight limit.
- 4-way stretch, vapour permeable, anti-bacterial welded PU cover with covered zip.
- Durable PVC base for a stable foundation.
- Suitable for very high risk with pump, high risk without pump.
- Fully modular system for ease of repair and cleaning.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HHMPRIMA.HPU250
Warranty	3 Years
Mattress Dimensions	2000 x 900 x 140 mm
Pump Dimensions	275 x 116 x 175 mm
Max. User Weight	247 kg / 39 st
Min. User Weight	N/A
Mattress Weight	9 kg
Pump Weight	2 kg
Number of Cells	10 alternating cells
Pressure Risk	Very High (with pump)
Cycle Time	1 in 2 over 10 mins
Medical Product	Class 1

PRODUCT CODES	
HHMPRIMA.HPU250	Prime Comfort Active Hybrid Mattress & 250 Pump

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: June 2026