



Prime Comfort Plus Static Mattress

The Prime Comfort Plus Static Mattress is an excellent static mattress providing great patient comfort. Side wall construction encases the foam, enhancing stability for the patient and assisting in patient transfer. A castellated foam top allows excess heat to escape from beneath the foam, preventing the patient from overheating whilst laid on the mattress. Keyhole slots in the high density base aid profiling, and it is available in three different sizes to suit a variety of patients.

KEY PRODUCT FEATURES

- ✓ **Side Wall Construction:** Encases the foam, enhancing stability and assisting with patient transfers.
- ✓ **Castellated Foam Top:** Allows excess heat to escape, preventing the patient from overheating.
- ✓ **Excellent Pressure Redistribution:** Suitable for patients at high risk of developing pressure ulcers.
- ✓ **Keyhole Slots in Base:** High density base with keyhole slots to aid profiling on electric beds.

HIGHLIGHTS

- ★ **Side Wall Construction**
- ★ **Castellated Foam Top**
- ★ **39 Stone Max Patient Weight**
- ★ **3 Sizes Available**
- ★ **3-Year Warranty**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- 4-way stretch, vapour permeable, water resistant PU cover with welded seams.
- Available in three sizes: 900mm, 1200mm and 1400mm widths. Designed for use on electric profiling bed frames.
- Side wall construction improves bed to chair transfers.
- Welded seams reduce fluid ingress and support infection control.
- Cover can be wiped down or machine washed up to 95°C.
- Also available as the Prime Comfort Active hybrid with 250 Pump.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HSMPRIMP
Warranty	3-Year Foam, 1-Year Cover
Mattress Dimensions	2000 x 900 x 140 mm
Max. User Weight	247 kg / 39 st
Pressure Risk	High
Foam Construction	Castellated top, high density base, side walls
Cover	4-way stretch PU, welded seams, Crib 5
Widths Available	900 mm / 1200 mm / 1400 mm
Profiling	Keyhole slots in base
Medical Product	Class 1

PRODUCT CODES	
HSMPRIMP	Reflect 2 1200 Static Mattress – 2000 x 1200 x 150 mm
HSMPRIMP.12	Prime Comfort Plus – 1200mm Wide
HSMPRIMP.14	Prime Comfort Plus – 1400mm Wide

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026