



Prime Comfort Static Mattress

The Prime Comfort Static Mattress is a quality mattress made of high resiliency foam, offering great pressure relief while remaining cost effective for any care environment. The foam is castellated for comfort and to facilitate air movement inside the mattress, reducing heat build-up, and the castellations are zonal for effective pressure relief. Suitable for patients at high risk of developing a pressure ulcer, up to 28 stone (177kg).

KEY PRODUCT FEATURES

- ✓ **Zonal Castellated Foam:** Castellations are zoned across the mattress for effective pressure relief.
- ✓ **High Resiliency Foam:** Quality foam construction giving great pressure relief at a cost effective price.
- ✓ **Improved Airflow:** Castellations facilitate air movement inside the mattress, reducing heat build-up.
- ✓ **Anti-Bacterial PU Cover:** Water resistant cover with welded seams and a Crib 5 fire rating.

HIGHLIGHTS

- ★ **Zonal Castellated Foam**
- ★ **High Risk Pressure Relief**
- ★ **28 Stone Max Patient Weight**
- ★ **Improved Airflow Reduces Heat Build-Up**
- ★ **3-Year Warranty**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Suitable for patients at high risk of developing a pressure ulcer.
Cost effective solution for any care environment.
- Designed for use on electric profiling bed frames.
- Welded seams reduce fluid ingress and support infection control.
- Additional sizes available, including 1200mm and 1400mm wide.
- Cover can be wiped down or machine washed up to 95°C.
- Also available as the Prime Comfort Active hybrid with 250 Pump.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HSMPRIM
Warranty	3-Year Foam, 1-Year Cover
Mattress Dimensions	2000 x 900 x 150 mm
Max. User Weight	177 kg / 28 st
Pressure Risk	High
Foam Construction	Zonal castellated high resiliency foam
Cover	4-way stretch PU, welded seams, Crib 5
Widths Available	900 mm / 1200 mm / 1400 mm
Mattress Weight	9 kg
Medical Product	Class 1

PRODUCT CODES	
HSMPRIM	Prime Comfort Static Mattress – 2000 x 900 x 150 mm
HSMPRIM.DV12	Prime Comfort Static Mattress – 1200mm Wide
HSMPRIM.14	Prime Comfort Static Mattress – 1400mm Wide

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026