



Reflect 2 1200 Static Mattress

The Reflect 2 1200 is a high grade foam profiling mattress designed for use on 1200mm wide beds, such as the Woburn Community 1200 and Ultimate 1200. A deep layer of castellated viscofoam rests upon a foam base; together the two foams evenly distribute patient weight and provide extra support across the additional width. The castellated top and keyhole slots aid profiling, and the no-turn design reduces the need to disturb the patient.

KEY PRODUCT FEATURES

- ✓ **Wide 1200mm Design:** Made for wider profiling beds including the Community 1200 and Ultimate 1200.
- ✓ **Castellated Viscofoam Top:** Memory foam effect evenly distributes weight and reduces pressure points.
- ✓ **Supportive Foam Base:** High resilience base prevents the mattress splitting and bottoming out.
- ✓ **Profiling Bed Compatible:** Castellated top and keyhole slots in the base aid profiling.

HIGHLIGHTS

- ★ **Fits 1200mm Wide Beds**
- ★ **High Risk Pressure Relief**
- ★ **39 Stone Max User Weight**
- ★ **No-Turn Design Low Maintenance**
- ★ **3-Year Warranty**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in [harvest-healthcare-ltd](#)
x [@HarvestHealthUK](#)

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- No monthly turning required, reducing patient disturbance.
- Castellated top allows excess heat to escape, preventing overheating.
- Suitable for prevention and treatment of pressure sores.
- 4-way stretch, vapour permeable, water resistant PU cover with welded seams.
- Crib 5 fire rated cover for care environment compliance. Cover can be wiped down or machine washed up to 95°C.
- Alternative treatment mattresses available from the Active range.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HFR147.12
Warranty	3-Year Foam, 1-Year Cover
Mattress Dimensions	2000 x 1200 x 150 mm
Max. User Weight	247 kg / 39 st
Pressure Risk	High
Foam Construction	Castellated viscofoam on foam base
Cover	4-way stretch PU, welded seams, Crib 5
Compatible Beds	Community 1200, Ultimate 1200
Turning Required	No
Medical Product	Class 1

PRODUCT CODES	
HFR147.12	Reflect 2 1200 Static Mattress – 2000 x 1200 x 150 mm

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026