



Reflect 2 Plus Static Mattress

The Reflect 2 Plus is Harvest Healthcare's bariatric static mattress for plus size patients up to 60 stone, in need of prevention and treatment for pressure sores. A deep layer of castellated viscoelastic foam and a firm foam base, combined with a sandwich layer of high load bearing foam, allow for greater prevention of pressure ulcers for bariatric patients. Available in 1200mm and 1400mm widths with 45 stone and 60 stone weight capacities.

KEY PRODUCT FEATURES

- ✓ **Bariatric Capacity:** 45 stone (285kg) and 60 stone (381kg) options for plus size patients.
- ✓ **High Load Bearing Foam:** Sandwich layer construction gives greater pressure ulcer prevention.
- ✓ **Castellated Viscofoam Top:** Evenly distributes weight, with keyhole slots to aid profiling.
- ✓ **Durable Bariatric Cover:** 4-way stretch PU with stitched seams, zipped to a PVC coated base.

HIGHLIGHTS

- ★ **Up to 60 Stone User Weight**
- ★ **1200mm & 1400mm Width Options**
- ★ **High Load Bearing Foam Sandwich**
- ★ **No-Turn Design Low Maintenance**
- ★ **3-Year Warranty**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in [harvest-healthcare-ltd](#)
x [@HarvestHealthUK](#)

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Silver-based antibacterial properties in the PU coating, washable and environmentally friendly.
- Navy blue cover; does not contain latex.
- No monthly turning required, reducing patient disturbance.
- Castellated top allows excess heat to escape, preventing overheating.
- Suitable for prevention and treatment of pressure sores. Cover can be wiped down or machine washed up to 95°C.
- Alternative bariatric treatment mattresses available from the Active range

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HFR227 (+ variants)
Warranty	3-Year Foam, 1-Year Cover
Widths Available	1200 mm / 1400 mm
Max. User Weight (45st versions)	285 kg / 45 st
Max. User Weight (60st versions)	381 kg / 60 st
Pressure Risk	High
Foam Construction	Viscofoam, high load bearing sandwich, foam base
Cover	4-way stretch PU, welded seams, Crib 5
Compatible Beds	Community 1200, Ultimate 1200
Turning Required	No
Medical Product	Class 1

PRODUCT CODES	
HFR147.12	Reflect 2 Plus – 1200mm, 45 Stone
HFR227.14	Reflect 2 Plus – 1400mm, 45 Stone
HFR227.60	Reflect 2 Plus – 1200mm, 60 Stone
HFR227.1400.60	Reflect 2 Plus – 1400mm, 60 Stone

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 4: July 2026