



Sandringham Active Mattress

The Sandringham Active Mattress and 250 Pump is a full mattress replacement system of air cells on a 50mm foam sub mattress. Made up of 18 individual cells, 16 active and 2 static pillow section cells, the foam sub mattress supports the patient in the event of power failure and adds stability. The compact digital pump provides Alternating, Constant Low Pressure, Max Inflate and Seated modes for patients at high risk, including existing tissue damage.

KEY PRODUCT FEATURES

- ✓ **18 TPU Cells:** 16 active and 2 static pillow section cells for stability and comfort.
- ✓ **50mm Foam Sub Mattress:** Supports the patient in the event of power failure and adds stability.
- ✓ **Full Replacement System:** Replaces the existing mattress for complete pressure care.
- ✓ **Versatile Pump Modes:** Alternating, Constant Low Pressure, Max Inflate and Seated modes.

HIGHLIGHTS

- ★ **50mm Foam Sub Mattress**
- ★ **High Risk Pressure Relief**
- ★ **Comfort Control Settings**
- ★ **200kg Max User Weight**

CONTACT US

☎ +44 (0)1709 377172
✉ sales@harvesthealthcare.co.uk

FOLLOW US

in harvest-healthcare-ltd
x @HarvestHealthUK

OUR LOCATION

📍 Sheaf House, Bradmarsh Business Park,
Bradmarsh Way, Rotherham S60 1BW

ADDITIONAL FEATURES

- Pillow section for patient stability and comfort.
- Max Inflate with auto return to alternating for patient safety.
- Seated mode to support the patient when the backrest is elevated.
- Pump maintains set pressures to reduce carer intervention.
- No minimum user weight limit.
- 4-way stretch, vapour permeable, anti-bacterial welded PU cover with 360° covered zip.
- Pump Auto Lock Out function to reduce risk of unintentional use.
- Fully modular system for ease of repair and cleaning.

For further information:

Scan the QR Code below to visit our website.



TECHNICAL DATA	
Product Code	HAMSAND001.HPU250
Warranty	3 Years
Mattress Dimensions	2000 x 900 x 160 mm
Pump Dimensions	275 x 116 x 175 mm
Max. User Weight	200 kg / 31 st
Min. User Weight	N/A
Total System Weight	11 kg
Number of Cells	18 (16 active, 2 static)
Pressure Risk	High
Cycle Time	1 in 2 over 10 mins
Medical Product	Class 1

PRODUCT CODES	
HAMSAND001.HPU250	Sandringham Active Mattress & 250 Pump

Important Safety Notice

Harvest Healthcare Ltd. supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour-permeable PU fabric that could cause a risk of suffocation in the event of prolonged facial contact. Mattress risk categories are to be used as a guide only. Risk assessments specific to the user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: June 2026